

Patient

Name: ROWE, BOYD P
Date of Birth: 02/06/1980
Sex: Male
Case Number: TN21-187332
Diagnosis: Adenocarcinoma, NOS

Specimen Information

Primary Tumor Site: Lung, NOS
Specimen Site: Liver
Specimen ID: C21-14187-A1
Specimen Collected: 02-Dec-2021
Test Report Date: 31-Dec-2021

Ordered By

Susanne Arnold, MD
UK Markey Cancer Center
 800 Rose Street, Whitney Hendrickson
 Building
 Lexington, KY 40536 USA
 (859) 257-4488

Results with Therapy Associations

BIOMARKER	METHOD	ANALYTE	RESULT	THERAPY ASSOCIATION		BIOMARKER LEVEL*
ALK	IHC	Protein	Positive 3+, 100%	BENEFIT	alectinib, ceritinib, crizotinib, lorlatinib	Level 1
	FISH	DNA-Tumor	Positive		brigatinib	Level 2
					alectinib, brigatinib, ceritinib, crizotinib, lorlatinib	Level 2
PD-L1 (22c3)	IHC	Protein	Positive, TPS: 85%	BENEFIT	cemiplimab, pembrolizumab	Level 1
PD-L1 (28-8)	IHC	Protein	Positive 1+, 75%	BENEFIT	nivolumab/ipilimumab combination	Level 1
BRAF	Seq	DNA-Tumor	Mutation Not Detected	LACK OF BENEFIT	dabrafenib and trametinib combination therapy, vemurafenib	Level 2
EGFR	Seq	DNA-Tumor	Mutation Not Detected	LACK OF BENEFIT	erlotinib, gefitinib	Level 2
KRAS	Seq	DNA-Tumor	Mutation Not Detected	LACK OF BENEFIT	sotorasib	Level 2
RET	Seq	RNA-Tumor	Fusion Not Detected	LACK OF BENEFIT	pralsetinib, selpercatinib	Level 2
ROS1	Seq	RNA-Tumor	Fusion Not Detected	LACK OF BENEFIT	entrectinib	Level 2

* Biomarker reporting classification: Level 1 – Companion diagnostic (CDx); Level 2 – Strong evidence of clinical significance or is endorsed by standard clinical guidelines; Level 3 – Potential clinical significance. Bolded benefit therapies, if present, highlight the most clinically significant findings.

The selection of any, all, or none of the matched therapies resides solely with the discretion of the treating physician. Decisions on patient care and treatment must be based on the independent medical judgment of the treating physician, taking into consideration all available information concerning the patient's condition, the FDA prescribing information for any therapeutic, and in accordance with the applicable standard of care. Whether or not a particular patient will benefit from a selected therapy is based on many factors and can vary significantly. All trademarks and registered trademarks are the property of their respective owners.

Important Note

An EML4-ALK variant 3b fusion was detected by RNA sequencing. This fusion is frequent in lung cancers (Bayliss 2016 Cell Mol Life Sci 73:1209). Intron 6 of EML4 (NM_019063.5) is joined to exon 20 of ALK (NM_004304.4) (Choi 2008 Cancer Res 68:4971). Different EML4-ALK variants may be associated with different clinical outcomes (Lin 2018 J Clin Oncol 36:1199).

The choice of ALK inhibitor treatment should be made with consideration of the patient's line of therapy and central nervous system (CNS) involvement. Alectinib, brigatinib, and lorlatinib are considered NCCN-preferred agents for the first-line setting. For patients with CNS involvement, alectinib, brigatinib, ceritinib, and lorlatinib have intracranial efficacy. Optimal sequencing of ALK-targeted therapy is an active area of investigation. Please see NCCN guidelines (NSCLC and CNS) and clinicaltrials.gov for more information.

This patient's tumor is positive for ALK and PD-L1 expression. For patients with driver alterations such as ALK, the use of ALK-targeted kinase inhibitors (TKIs) should be prioritized. Clinical trials are exploring whether there is a role for immune checkpoint inhibitors after progression on available ALK TKIs. Please see NCCN guidelines and clinicaltrials.gov for more information.

Please note that multiple companion diagnostic assays (antibodies) have been utilized to assess PD-L1 expression. Each test has different performance characteristics, therefore, the results will not always be concordant.

Cancer-Type Relevant Biomarkers

Biomarker	Method	Analyte	Result
ALK	Seq	RNA-Tumor	Pathogenic Fusion
		DNA-Tumor	Mutation Not Detected
MSI	Seq	DNA-Tumor	Stable
Mismatch Repair Status	IHC	Protein	Proficient
NTRK1/2/3	Seq	RNA-Tumor	Fusion Not Detected
Tumor Mutational Burden	Seq	DNA-Tumor	Low, 2 mut/Mb
BRAF	Seq	RNA-Tumor	Fusion Not Detected
ERBB2 (Her2/Neu)	Seq	DNA-Tumor	Mutation Not Detected
FGFR3	Seq	RNA-Tumor	Fusion Not Detected
KEAP1	Seq	DNA-Tumor	Mutation Not Detected
KRAS	CNA-Seq	DNA-Tumor	Amplification Not Detected

Biomarker	Method	Analyte	Result
MET	Seq	RNA-Tumor	Variant Transcript Not Detected
	CNA-Seq	DNA-Tumor	Amplification Not Detected
	Seq	DNA-Tumor	Mutation Not Detected
NFE2L2	Seq	DNA-Tumor	Mutation Not Detected
NRG1	Seq	RNA-Tumor	Fusion Not Detected
PD-L1 (SP142)	IHC	Protein	Negative, IC: 0% Negative, TC: 0%
RB1	Seq	DNA-Tumor	Mutation Not Detected
RET	Seq	DNA-Tumor	Mutation Not Detected
STK11	Seq	DNA-Tumor	Mutation Not Detected
TP53	Seq	DNA-Tumor	Mutation Not Detected

Genomic Signatures

Biomarker	Method	Analyte	Result
Microsatellite Instability (MSI)	Seq	DNA-Tumor	Stable
Tumor Mutational Burden (TMB)	Seq	DNA-Tumor	Result: Low 2 Low10High
Genomic Loss of Heterozygosity (LOH)	Seq	DNA-Tumor	Low - 7% of tested genomic segments exhibited LOH (assay threshold is $\geq 16\%$)

Genes Tested with Pathogenic or Likely Pathogenic Alterations

Gene	Method	Analyte	Variant Interpretation	Protein Alteration	Exon	DNA Alteration	Variant Frequency %
ALK	Seq	RNA-Tumor	Pathogenic Fusion	EML4-ALK	20	-	-
CDKN2A	Seq	DNA-Tumor	Pathogenic Variant	p.R80*	2	c.238C>T	24

Unclassified alterations for DNA and RNA sequencing can be found in the MI Portal.
Formal nucleotide nomenclature and gene reference sequences can be found in the Appendix of this report.

Human Leukocyte Antigen (HLA) Genotype Results

The impact of HLA genotypes on drug response and prognosis is an active area of research. These results can help direct patients to clinical trials recruiting for specific genotypes. Please see www.clinicaltrials.gov for more information.

Gene	Method	Analyte	Genotype
MHC CLASS I			
HLA-A	Seq	RNA-Tumor	A*29:02,A*33:01
HLA-B	Seq	RNA-Tumor	B*44:04,B*14:02
HLA-C	Seq	RNA-Tumor	C*16:01,C*08:02

Please note that the HLA sequencing data above was obtained from expressed RNA transcripts and not from DNA sequencing reads.
HLA sequencing data from DNA may be different from what is reported here.
HLA genotypes with only one allele are either homozygous or have loss-of-heterozygosity at that position.

Immunohistochemistry Results

Biomarker	Result	Biomarker	Result
ALK	Positive 3+, 100%	PD-L1 (28-8)	Positive 1+, 75%
MLH1	Positive 2+, 85%	PD-L1 (SP142)	Negative, IC: 0% Negative, TC: 0%
MSH2	Positive 2+, 90%	PMS2	Positive 2+, 85%
MSH6	Positive 2+, 85%	PTEN	Positive 1+, 100%
PD-L1 (22c3)	Positive, TPS: 85%		

In Situ Hybridization Results

Fusion Detected
ALK

Genes Tested with Indeterminate Results by Tumor DNA Sequencing

COL2A1	NOTCH1	PIK3R2	PTPN11	PTPRD	RAC1	RASA1	STAG2	TERT	TRAF7	XRCC1	XRCC2
MED12	PIK3CB										

Genes in this table were ruled indeterminate due to low coverage for some or all exons.

The results in this report were curated to represent biomarkers most relevant for the submitted cancer type. These include results important for therapeutic decision-making, as well as notable alterations in other biomarkers known to be involved in oncogenesis. Additional results, including genes with normal findings and unclassified alterations can be found in the MI Portal at miportal.carismolecularintelligence.com. If you do not have an MI Portal account, or need assistance accessing it, please contact Caris Customer Support at (888) 979-8669.

Notes of Significance

SEE APPENDIX FOR DETAILS

An EML4-ALK variant 3b fusion was detected by RNA sequencing. This fusion is frequent in lung cancers (Bayliss 2016 Cell Mol Life Sci 73:1209). Intron 6 of EML4 (NM_019063.5) is joined to exon 20 of ALK (NM_004304.4) (Choi 2008 Cancer Res 68:4971). Different EML4-ALK variants may be associated with different clinical outcomes (Lin 2018 J Clin Oncol 36:1199). 12/22/2021

Clinical Trials Connector™ opportunities based on biomarker expression: 144 Targeted Therapy Trials. See page 6 for details.

Specimen Information

Specimen ID: C21-14187-A1

Specimen Collected: 12/02/2021

Specimen Received: 12/15/2021

Testing Initiated: 12/15/2021

Gross Description: 1 (A) Paraffin Block - Client ID(C21-14187-A1) from UK Healthcare Albert B Chandler Hospital, Lexington, KY, with the corresponding pathology report labeled "C21-14187".

Pathologic Diagnosis: Liver, right lobe, CT-guided core biopsy with touch preps: Positive for malignancy, adenocarcinoma compatible with lung primary.

Dissection Information: Molecular testing of this specimen was performed after harvesting of targeted tissues with an approved manual microdissection technique. Candidate slides were examined under a microscope and areas containing tumor cells (and separately normal cells, when necessary for testing) were circled. A laboratory technician harvested targeted tissues for extraction from the marked areas using a dissection microscope.

Clinical Trials Connector™

For a complete list of open, enrolling clinical trials visit MI Portal to access the [Clinical Trials Connector](#). This personalized, real-time web-based service provides additional clinical trial information and enhanced searching capabilities, including, but not limited to:

- Location: filter by geographic area
- Biomarker(s): identify specific biomarkers associated with open clinical trials to choose from
- Drug(s): search for specific therapies
- Trial Sponsor: locate trials based on the organization supporting the trial(s)

The Clinical Trials Connector lists agents that are matched to available clinical trials according to biomarker status. In some instances, older-generation agents may still be relevant in the context of new combination strategies and, therefore, will still appear on this report.

Visit www.Carismolecularintelligence.com to view all matched trials. Therapeutic agents listed below may or may not be currently FDA approved for the tumor type tested.

TARGETED THERAPY CLINICAL TRIALS (144)				
Drug Class	Biomarker	Method	Analyte	Investigational Agent(s)
ALK inhibitors (23)	ALK	IHC	Protein	alectinib, brigatinib
	ALK	RNA-Seq	RNA-Tumor	
Immunomodulatory agents (121)	PD-L1	IHC	Protein	CDX-527, atezolizumab, avelumab, camrelizumab, cemiplimab

() = represents the total number of clinical trials identified by the Clinical Trials Connector for the provided drug class or table.

The Clinical Trials Connector may include trials that enroll patients with additional screening of molecular alterations. In some instances, only specific gene variants may be eligible.

Disclaimer

Decisions on patient care and treatment must be based on the independent medical judgment of the treating physician, taking into consideration all available information concerning the patient's condition, the FDA prescribing information for any therapeutic, and in accordance with the applicable standard of care. Drug associations provided in this report do not guarantee that any particular agent will be effective for the treatment of any patient or for any particular condition. Caris Life Sciences expressly disclaims and makes no representation or warranty whatsoever relating, directly or indirectly, to the performance of services, including any information provided and/or conclusions drawn from therapies that are included or omitted from this report. Whether or not a particular patient will benefit from a selected therapy is based on many factors and can vary significantly. The selection of therapy, if any, resides solely in the discretion of the treating physician.

Individual assays that are available through Caris Molecular Intelligence® include both Laboratory Developed Tests (LDT) and U.S. Food and Drug Administration (FDA) approved or cleared tests. Caris MPI, Inc. d/b/a Caris Life Sciences® is certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical laboratory testing, including all of the assays that comprise the Caris Molecular Intelligence®. The LDTs were developed and their performance characteristics determined by Caris. The LDTs have not been cleared or approved by the U.S. Food and Drug Administration. Caris' CLIA certification number is located at the bottom of each page of this report. Certain tests have not been cleared or approved by the FDA. The FDA has determined that clearance or approval is not necessary for certain laboratory developed tests. Caris LDTs are used for clinical purposes. They are not investigational or for research.

The information presented in the Clinical Trials Connector™ section of this report, if applicable, is compiled from sources believed to be reliable and current. However, the accuracy and completeness of the information provided herein cannot be guaranteed. The clinical trials information present in the biomarker description was compiled from www.clinicaltrials.gov. The contents are to be used only as a guide, and health care providers should employ their best comprehensive judgment in interpreting this information for a particular patient. Specific eligibility criteria for each clinical trial should be reviewed as additional inclusion criteria may apply.

All materials, documents, data, data software, information and/or inventions supplied to customers by or on behalf of Caris or created by either party relating to the services shall be and remain the sole and exclusive property of Caris. Customer shall not use or disclose the information provided by Caris through the services or related reports except in connection with the treatment of the patient for whom the services were ordered and shall not use such property for, or disseminate such property to, any third parties without expressed written consent from Caris. Customer shall deliver all such property to Caris immediately upon demand or upon Caris ceasing to provide the services.

Caris Molecular Intelligence is subject to Caris' intellectual property. Patent www.carislifesciences.com/ip.

A handwritten electronic signature in black ink, appearing to read "Mylinh Mac".

Electronic Signature

Mylinh Mac, MD

12/31/2021

Mutational Analysis by Next-Generation Sequencing (NGS)

TUMOR MUTATIONAL BURDEN	
Mutations / Megabase	Result
2	Low

TMB Methods

Tumor Mutational Burden (TMB) analysis was performed based on Next Generation Sequencing analysis from genomic DNA isolated from a formalin-fixed paraffin embedded tumor sample using the Illumina platform. TMB is calculated using nonsynonymous, in-frame indel, and frameshift indel mutations that have not been previously reported as germline alterations in the Genome Aggregation Database (gnomAD) and dbSNP151 or as common benign variants identified by Caris geneticists. The cutoff of 10 mutations/megabase was established in NSCLC and its applicability towards other tumor types has not been established at this time. Preliminary data supporting the pan-tumor FDA approval are available in Marabelle et al., "Association of tumour mutational burden in patients with advanced solid tumours treated with pembrolizumab: prospective biomarker analysis of the multicohort, open-label, phase 2 KEYNOTE-158 study". Lancet Oncology, 2020. Caris Life Sciences is a participant in the Friends of Cancer Research TMB Harmonization Project. (Merino et al., 2020).

MICROSATELLITE INSTABILITY ANALYSIS		
Test	Interpretation	Result
MSI	No microsatellite instability detected.	Stable
	Procedure: NGS	

Microsatellite Instability Analysis

Microsatellite instability status by NGS (MSI-NGS) is measured by the direct analysis of known microsatellite regions sequenced in the CMI NGS panel. To establish clinical thresholds, MSI-NGS results were compared with results from over 2,000 matching clinical cases analyzed with traditional, PCR-based methods. Genomic variants in the microsatellite loci are detected using the same depth and frequency criteria as used for mutation detection. Only insertions and deletions resulting in a change in the number of tandem repeats are considered in this assay. Some microsatellite regions with known polymorphisms or technical sequencing issues are excluded from the analysis. The total number of microsatellite alterations in each sample are counted and grouped into three categories: High (≥ 116 MSI loci altered), Equivocal (113-115 MSI loci altered) and Stable (≤ 112 MSI loci altered).

GENOMIC LOSS OF HETEROZYGOSITY	
Test	Result
Genomic Loss of Heterozygosity (LOH)	Low - 7% of tested genomic segments exhibited LOH (assay threshold is $\geq 16\%$)

Genomic Loss of Heterozygosity Analysis:

In order to calculate genomic loss-of-heterozygosity (LOH), the 22 autosomal chromosomes are split into 552 segments and the LOH of single nucleotide polymorphisms (SNPs) within each segment is calculated. Caris WES data consist of approximately 250k SNPs spread across the genome. SNP alleles with frequencies skewed towards 0 or 100% indicate LOH (heterozygous SNP alleles have a frequency of 50%). In this assay, a segment is determined to have LOH if the average SNP variant frequency is skewed more than $\pm 15\%$ from the heterozygous frequency of 50% (p-value < 0.02 after correction vs. a negative control). The final call of genomic LOH is based on the percentage of all 552 segments with observed LOH (High $\geq 16\%$, Low $< 16\%$; if fewer than 3,000 SNPs can be read, the test is reported as Indeterminate). A normal epithelial ovarian genome (NA12878) that has no non-polymorphic variants, gene fusions or other cancer hallmarks, is used as a negative control. Segment sizes range from 2-6 Mb, depending on segment proximity to the centromeres or telomeres. 99% of segments are at least 5Mb. Segments excluded from the calculation of genomic LOH include those spanning $\geq 90\%$ of a whole chromosome or chromosome arm and segments which are not covered by the SNP backbone and the WES panel. The 250k SNPs consist of 200K from exonic regions and 50K from intronic regions, with a minimum of 17 SNPs per Mb of genome sequence.

Additional Next-Generation Sequencing results continued on the next page. >

PATIENT: ROWE, BOYD P (06-FEB-1980)

TN21-187332

PHYSICIAN: SUSANNE ARNOLD, MD

Mutational Analysis by Next-Generation Sequencing (NGS)

Gene	Analyte	Variant Interpretation	Protein Alteration	Exon	DNA Alteration	Variant Frequency %	Transcript ID
CDKN2A	DNA-Tumor	Pathogenic Variant	p.R80*	2	c.238C>T	24	NM_000077.4

Interpretation: A pathogenic mutation was detected in CDKN2A (p16).

CDKN2A, or cyclin-dependent kinase inhibitor 2A, is a tumor suppressor gene that encodes two cell-cycle regulatory proteins, p16INK4A and p14ARF. As upstream regulators of the retinoblastoma (RB) and p53 signaling pathways, CDKN2A controls the induction of cell cycle arrest in damaged cells that allows for repair of DNA. Loss of CDKN2A through whole-gene deletion, point mutation, or promoter methylation leads to disruption of these regulatory proteins and consequently dysregulation of growth control. Somatic CDKN2A mutations are documented to occur in squamous cell lung cancers, head and neck cancer, colorectal cancer, chronic myelogenous leukemia and malignant pleural mesothelioma. Currently, there are agents that target downstream of CDKN2A such as CDK4/6 inhibitors which function by restoring the cell's ability to induce cell cycle arrest. Germline CDKN2A mutations are associated with melanoma-pancreatic carcinoma syndrome, which increases the risk for familial malignant melanoma and pancreatic cancer.

Additional Next-Generation Sequencing results continued on the next page. >

PATIENT: ROWE, BOYD P (06-FEB-1980)

TN21-187332

PHYSICIAN: SUSANNE ARNOLD, MD

Mutational Analysis by Next-Generation Sequencing (NGS)

GENES TESTED WITH INDETERMINATE * RESULTS BY TUMOR DNA SEQUENCING

COL2A1	PIK3CB	PTPRD	STAG2	XRCC1	
MED12	PIK3R2	RAC1	TERT	XRCC2	
NOTCH1	PTPN11	RASA1	TRAF7		

* Genes in this table were ruled indeterminate due to low coverage for some or all exons.

For a complete list of genes tested, visit www.Carismolecularintelligence.com/profilemenu.

NGS Methods

Next Generation Sequencing for WES (Whole Exome Sequencing): Direct sequence analysis was performed on genomic DNA isolated from a micro-dissected, formalin-fixed paraffin-embedded tumor sample using the Illumina NovaSeq 6000 sequencers. A hybrid pull-down panel of baits designed to enrich for more than 700 clinically relevant genes at high coverage and high read-depth was used, along with another panel designed to enrich for an additional >20,000 genes at lower depth. A 500Mb SNP backbone panel (Agilent Technologies) was added to assist with gene amplification/deletion measurements and other analyses. The performance of the CMI WES assay was validated for sequencing variants, copy number alteration, tumor mutational burden and micro-satellite instability. The test was validated to 50ng of input and has a PPV of 0.99 against a previously validated NGS assay. CMI WES can detect variants with tumor nuclei as low as 20%, and will detect variants down to 5% variant frequency with an average depth of at least 500x. This test has a sensitivity to detect as low as approximately 10% population of cells containing a mutation in all exons from the high read-depth clinical genes and 99% of all exons in the 20K whole exome regions. CMI WES is currently validated to detect <44bp indels. The reference genome for the transcript ID is hg38 with hg19 liftOver calculations performed for the high read-depth gene panel. While the vast majority of exons in the exome are covered by the assay, technical constraints preclude the coverage of every exon. Of the high read-depth genes with the most relevance to cancer, the following have only partial exon coverage: ARID1B, ASXL2, CDH23, CDKN1C, CHEK2, CYP2D6, DIS3L2, EIF1AX, FAT3, FLT4, FOXO3, HSP90AA1, HSP90AB1, KMT2C, MAGI2, MAML2, MDS2, MLLT3, NCOR1, NOTCH2, NSD3, PDE4DIP, PMS2, RAC1, RAD52, RANBP2, RB1, RHEB, RPL10, RPL22, SBDS, SET, SMC3, SRSF3, STAT5B, SUZ12, TCEA1, TOP3B, TSHZ3, USP6, and ZFXH3. For a complete list of what is covered, please contact Caris Customer Support.

Copy Number Alterations by Next-Generation Sequencing (NGS)

CNA Methods

The copy number alteration (CNA) of each exon is determined by a calculation using the average sequencing depth of the sample along with the sequencing depth of each exon and comparing this calculated result to a pre-calibrated value. If all exons within the gene of interest have an average of ≥ 3 copies and the average copy number of the entire gene is ≥ 6 copies, the gene result is reported as amplified. If an average of ≥ 4 but < 6 copies of a gene are detected, or if the average copy number of the gene is ≥ 6 copies but contains exons with an average of < 3 copies, the gene result is reported as intermediate. If a gene has > 0.1 copies but < 1.3 copies, the gene result is reported as deleted. If an average > 1.9 copies but < 4 copies of a gene are detected, the gene result is reported as no amplification detected or deletion not detected. A complete list of copy number alteration genes are available upon request.

Gene Fusion and Transcript Variant Detection by RNA Sequencing

GENES TESTED WITH GENE FUSION OR TRANSCRIPT VARIANT DETECTED

Biomarker	Fusion/Isoform	Splice Site	Transcript ID	Variant Interpretation
ALK	EML4:ALK	intron 6:exon 20	NM_019063/NM_004304	Pathogenic Fusion

Interpretation: An EML4-ALK variant 3b fusion was detected by RNA sequencing. This fusion is frequent in lung cancers (Bayliss 2016 Cell Mol Life Sci 73:1209). Intron 6 of EML4 (NM_019063.5) is joined to exon 20 of ALK (NM_004304.4) (Choi 2008 Cancer Res 68:4971). Different EML4-ALK variants may be associated with different clinical outcomes (Lin 2018 J Clin Oncol 36:1199).

Whole Transcriptome Sequencing (WTS) Methods

Gene fusion and variant transcript detection, including HLA genotyping, were performed on mRNA isolated from a formalin-fixed paraffin-embedded tumor sample using the Agilent SureSelectXT Low Input Library prep chemistry, optimized for FFPE tissue, in conjunction with the SureSelect Human All Exon V7 bait panel (48.2 Mb) and the Illumina NovaSeq. In addition to transcript variants, this assay is designed to detect fusions occurring at known and novel breakpoints within genes. Only a portion of genes tested are included in this report. The genes included in this report represent the subset of genes most commonly associated with cancer. All results can be provided by request. For fusions and non-HLA variant transcripts, analytical validation of this test demonstrated $\geq 97\%$ Positive Percent Agreement (PPA), $\geq 99\%$ Negative Percent Agreement (NPA) and $\geq 99\%$ Overall Percent Agreement (OPA) with a validated comparator method. For HLA genotyping, analytical validation of this test demonstrated $\geq 99\%$ Positive Percent Agreement (PPA), $\geq 98\%$ Negative Percent Agreement (NPA) and $\geq 99\%$ Overall Percent Agreement (OPA) with a validated comparator method. The versioned reference identifier used for the transcript ID was Feb.2009 (GRCh37/hg19). The complete list of unclassified alterations for RNA Whole Transcriptome Sequencing are available by request. HLA results are not available in New York State.

Protein Expression by Immunohistochemistry (IHC)

Biomarker	Patient Tumor			Thresholds
	Staining Intensity (0, 1+, 2+, 3+)	Percent of cells	Result	Conditions for a Positive Result:
MLH1	2 +	85	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained
MSH2	2 +	90	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained
MSH6	2 +	85	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained
PMS2	2 +	85	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained
PTEN	1 +	100	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained

PD-L1 TUMOR CELL STAINING				
Biomarker	Patient Tumor			Thresholds
	Staining Intensity (0, 1+, 2+, 3+)	Percent of cells	Result	Conditions for a Positive Result:
PD-L1 (28-8)	1 +	75%	Positive	Intensity $\geq 1+$ and $\geq 1\%$ of cells stained
PD-L1 (SP142)	0	100%	Negative	Intensity $\geq 1+$ and $\geq 50\%$ of cells stained

PD-L1 (28-8): Scoring was based on percentage of viable tumor cells showing partial or complete membrane staining at any intensity.

PD-L1 (SP142): TC scoring was based on the presence of discernible PD-L1 membrane staining of any intensity in $\geq 50\%$ of viable tumor cells.


PD-L1 TUMOR PROPORTION SCORE (TPS)			
Biomarker	Result	TPS	Threshold
PD-L1 (22c3)	Positive	85%	TPS $\geq 1\%$

PD-L1 22c3: Scoring was based on the percentage of viable tumor cells showing partial or complete membrane staining. There are three categories of PD-L1 expression defined by the PD-L1 22c3 IHC pharmDx NSCLC interpretation guide: TPS $< 1\%$ (negative), TPS $\geq 1\%$ and TPS $\geq 50\%$.

PD-L1 IMMUNE CELL (IC) SCORE			
Biomarker	Result	IC	Threshold
PD-L1 (SP142)	Negative	0%	$\geq 10\%$

PD-L1 (SP142): IC scoring was based on discernible PD-L1 staining of any intensity in tumor-infiltrating immune cells covering $\geq 10\%$ of tumor area occupied by tumor cells, associated intratumoral or contiguous peritumoral stroma.

Clones used: MLH1 (M1), MSH2 (G219-1129), MSH6 (SP93), PMS2 (A16-4), PD-L1 (SP142), PD-L1 (22c3), PTEN (6H2.1), PD-L1 (28-8).



Electronic Signature
James Lee, MD, PhD
12/19/2021

Additional IHC results continued on the next page. >

PATIENT: ROWE, BOYD P (06-FEB-1980)

TN21-187332

PHYSICIAN: SUSANNE ARNOLD, MD

Protein Expression by Immunohistochemistry (IHC)

Biomarker	Patient Tumor			Thresholds
	Staining Intensity (0, 1+, 2+, 3+)	Percent of cells	Result	Conditions for a Positive Result:
ALK	3 +	100	Positive	Intensity $\geq 3+$ and $\geq 1\%$ of cells stained

Clones used: ALK (D5F3).



Electronic Signature
Mylinh Mac, MD
12/31/2021

IHC Methods

The Laboratory Developed Tests (LDT) immunohistochemistry (IHC) assays were developed and their performance characteristics determined by Caris Life Sciences. These tests have not been cleared or approved by the US Food and Drug Administration. The FDA has determined that such clearance or approval is not currently necessary. Interpretations of all immunohistochemistry (IHC) assays were performed manually by a board certified pathologist using a microscope and/or digital whole slide image(s).

The following IHC assays were performed using FDA-approved companion diagnostic or FDA-cleared tests consistent with the manufacturer's instructions: ALK (VENTANA ALK (D5F3) CDx Assay, Ventana), ER (CONFIRM anti-Estrogen Receptor (ER) (SP1), Ventana), PR (CONFIRM anti-Progesterone Receptor (PR) (1E2), Ventana), HER2/neu (PATHWAY anti-HER-2/neu (4B5), Ventana), Ventana, PD-L1 22c3 (pharmDx, Dako), PD-L1 SP142 (VENTANA, Ventana in urothelial carcinomas, breast carcinoma and non-small cell lung cancer; drug association only in urothelial, triple negative breast cancers and non-small cell lung cancer), and PD-L1 28-8 (pharmDx, Dako).

HER2 results and interpretation follow the ASCO/CAP scoring criteria.

Rearrangement by fluorescence in situ Hybridization (FISH)

Probe Name	Result	ISCN/Probe Description
ALK	Fusion Detected	ALK (2p23.2) Break Apart
<i>Reference Range:</i> Positivity for ALK rearrangement is defined as >25 positive cells out of the 50 cells analyzed. A sample is considered negative if <5 positive cells are present out of the 50 cells analyzed. In cases where 5-25 cells are positive, the sample is considered equivocal, and an additional 50 cells are analyzed by a second technologist. From this expanded analysis, if ≥15 cells out of the 100 cells analyzed are positive for ALK rearrangement, the sample is considered positive. If <15 positive cells are observed out of the 100 analyzed, the sample is considered negative.		



Electronic Signature

Mylinh Mac, MD

12/31/2021

FISH Methods

The reflex process that Phenopath Laboratories (Seattle, WA) utilizes has a different reference range than that listed above. Please refer to the transcribed methodology below. CMPI will retain a copy of the original report methodology for this case.

FISH was performed on formalin fixed paraffin embedded tissue sections using an ALK Break Apart Probe for the ALK IVD (2p23)(Abbott Molecular, Abbott Park, III). A total of 50 cells were scored 16.0% of cells had a signal pattern consistent with a rearrangement involving the ALK gene at 2p23.

Each ALK test performed is evaluated by a pathologist including analysis of an H&E-stained section for tumor cell adequacy. Probed slides are scanned in areas selected by a pathologist using the MetaSystems™ imaging platform. Tumor cells are then scored in two steps. In the first step using an ALK 2p23 probe set (Biocare Medical, Concord, CA), denoted in the table as ALK (2p23), 50 tumor cells are evaluated. If < 10% of the tumor cells exhibit a breakapart or single red pattern, the case is negative. If ≥10% exhibit a breakapart or single red pattern, the case is reflexed for testing with the Vysis® ALK IVD Breakapart kit (Abbott Molecular, Abbott Park, III), denoted in the table as ALK IVD (2p23), using the ALK IVD scoring protocol. The threshold for positivity is 15% and only the Vysis® ALK IVD result is reported when performed. In select cases only the Vysis® ALK IVD Breakapart kit is used. The performance characteristics of this test have been validated by PhenoPath Laboratories.

ALK 2p23 was developed and its performance characteristics determined by PhenoPath Laboratories. This laboratory is certified under the Clinical Laboratory Improvement Amendments of 1988 as qualified to perform high complexity clinical testing. Testing was performed by PhenoPath Laboratories, 551 North 34th Street, Suite 100, Seattle, WA 98103-8675. The specimen was tested on 21 DEC 2021 and results were reported out on 30 DEC 2021. This test was interpreted by B. Jennifer Kum, M.D., Pathologist.

References

#	Drug	Biomarker	Reference
1	alectinib	ALK	Camidge, D.R., A.T. Shaw, et al. (2019). "Updated Efficacy and Safety Data and Impact of the EML4-ALK Fusion Variant on the Efficacy of Alectinib in Untreated ALK-Positive Advanced Non-Small Cell Lung Cancer in the Global Phase III ALEX Study." J Thorac Oncol 14(7): 1233-1243. View Citation Online
2	alectinib	ALK	Gadgeel, S., D.R. Camidge, et al. (2018). "Alectinib versus crizotinib in treatment-naïve anaplastic lymphoma kinase-positive (ALK+) non-small-cell lung cancer: CNS efficacy results from the ALEX study." Ann Oncol 29 (11): 2214-2222. View Citation Online
3	alectinib, brigatinib, ceritinib, crizotinib, lorlatinib	ALK	National Comprehensive Cancer Network. NCCN Clinical Practice Guidelines in Oncology. Non-Small Cell Lung Cancer Version 1.2020
4	alectinib, brigatinib, ceritinib, crizotinib, lorlatinib	ALK	Lindeman, N.I., Y. Yatabe, et al. (2018). "Updated Molecular Testing Guideline for the Selection of Lung Cancer Patients for Treatment With Targeted Tyrosine Kinase Inhibitors Guideline From the College of American Pathologists, the International Association for the Study of Lung Cancer, and the Association for Molecular Pathology." J Thorac Oncol 13(3): 323-358. View Citation Online
5	entrectinib	ROS1	Desai AV, Brodeur GM, Foster J, et al. Phase I study of entrectinib (RXDX-101), a TRK, ROS1, and ALK inhibitor, in children, adolescents, and young adults with recurrent or refractory solid tumors. J Clin Oncol. 2018;36 (suppl;abstr 10536). doi: 10.1200/JCO.2018.36.15_suppl.10536. View Citation Online
6	entrectinib	ROS1	Demetri, G.D., R.D., Doebele, et al. (2018). "Efficacy and safety of entrectinib in patients with NTRK fusion-positive tumors: pooled analysis of STARTRK-2, STARTRK-1 and ALKA-372-001. Presented at: 2018 ESMO Congress; October 19-23, 2018; Munich, Germany. Abstract LBA17. View Citation Online
7	crizotinib	ALK	Solomon, B. J., T. S. Mok, et al. (2018). "Final overall survival analysis from a study comparing first-line crizotinib versus chemotherapy in ALK-mutation positive Non-small-cell-lung cancer." J Clin Oncol 36: 2251-2258. View Citation Online
8	crizotinib	ALK	van der Wekken, A. J., H.J.M. Groen, et al. (2017). "Dichotomous ALK IHC is a better predictor for ALK inhibition outcome than traditional ALK FISH in advanced Non-small cell lung cancer." Clin Cancer Res 23(15): 4251-4258. View Citation Online
9	brigatinib, crizotinib	ALK	Thorne-Nuzzo, T., P. Towne, et al. (2017). "A Sensitive ALK Immunohistochemistry Companion Diagnostic Test Identifies Patients Eligible for Treatment with Crizotinib." J Thorac Oncol 12(5): 804-813. View Citation Online
10	brigatinib	ALK	Camidge, D.R., S. Popat, et al. (2018). "Brigatinib versus Crizotinib in ALK-Positive Non-Small-Cell Lung Cancer." N Engl J Med 379(21): 2027-2039. View Citation Online
11	brigatinib	ALK	Lin, J.L., G.J. Riely, et al. (2018). "Brigatinib in Patients with Alectinib-refractory ALK-positive NSCLC." J Thorac Onc 13(10): 1530-1538. View Citation Online
12	brigatinib	ALK	Reckamp, K., J. Lee, et al. (2019). "Comparative efficacy of brigatinib versus ceritinib and alectinib in patients with crizotinib-refractory anaplastic lymphoma kinase-positive non-small cell lung cancer." Curr Med Res Opin 35(4):569-576. View Citation Online
13	ceritinib, lorlatinib	ALK	Soria, J.C., de Castro G., et al. (2017). "First-line ceritinib versus platinum-based chemotherapy in advanced ALK-rearranged non-small-cell lung cancer (ASCEND-4): a randomised, open-label, phase 3 study." Lancet 389: 917-929. View Citation Online
14	ceritinib, lorlatinib	ALK	Solomon, B. J., A. T. Shaw, et al. (2018). "Lorlatinib in patients with ALK-positive non-small cell lung cancer: results from a global phase 2 study." Lancet Oncol 19:1654-1667. View Citation Online
15	ceritinib, lorlatinib	ALK	Shaw, A.T., E. Felip, et al. (2017). "Ceritinib versus chemotherapy in patients with ALK-rearranged non-small-cell lung cancer previously given chemotherapy and crizotinib (ASCEND-5): a randomised, controlled, open-label, phase 3 trial." Lancet Oncol 18 (7):874-886. View Citation Online
16	nivolumab/ipilimumab combination	PD-L1 (28-8)	Hellmann, M.D., S.S. Ramalingam, et al., (2019). "Nivolumab plus Ipilimumab in Advanced Non-Small-Cell Lung Cancer." N Engl J Med 381:2020-31. View Citation Online

References

#	Drug	Biomarker	Reference
17	pralsetinib, selpercatinib	RET	National Comprehensive Cancer Network. NCCN Clinical Practice Guidelines in Oncology. Non-Small Cell Lung Cancer Version 6.2020 View Citation Online
18	selpercatinib	RET	Drilon, A., V. Subbiah, et al. (2020). "Efficacy of Selpercanib in RET Fusion-Positive Non-Small-Cell Lung Cancer." N Engl J Med. 383 (9): 813-824. View Citation Online
19	pralsetinib	RET	Gainor J.F., V. Subbiah, et al. (2020). "Registrational dataset from the phase I/II ARROW trial of pralsetinib (BLU-667) in patients (pts) with advanced RET fusion+ non-small cell lung cancer (NSCLC)." J Clin Oncol. 38(suppl):9515. View Citation Online
20	cemiplimab	PD-L1 (22c3)	Sezer, A., P. Rietschel, et al., (2020). "EMPOWER-Lung 1: Phase III first-line (1L) cemiplimab monotherapy vs platinum-doublet chemotherapy (chemo) in advanced non-small cell lung cancer (NSCLC) with programmed cell death-ligand 1 (PD-L1) ≥50%." Ann Oncol 31 (suppl_4): S1142-S1215 View Citation Online
21	sotorasib	KRAS	Hong, D.S., B.T. Li, et al. (2020). "KRAS G12C Inhibition with Sotorasib in Advanced Solid Tumors." N Engl J Med. 383(13): 1207-1217. View Citation Online
22	sotorasib	KRAS	Li, B.T., J. Wolf, et al. (2021). "CodeBreak 100: Registrational Phase 2 Trial of Sotorasib in KRAS p.G12C Mutated Non-small Cell Lung Cancer". J Thorac Oncol. 16 (3): S61. View Citation Online
23	erlotinib, gefitinib	EGFR	Maemondo, M., T. Nukiwa, et al. (2010). "Gefitinib or chemotherapy for non-small cell lung cancer with mutated EGFR." N. Engl. J. Med. 362:2380-8. View Citation Online
24	erlotinib, gefitinib	EGFR	Brugger, W., F. Cappuzzo, et al. (2011). "Prospective molecular marker analyses of EGFR and KRAS from a randomized, placebo-controlled study of erlotinib maintenance therapy in advanced non-small-cell lung cancer." J. Clin. Oncol. 29:4113-4120. View Citation Online
25	erlotinib, gefitinib	EGFR	Keedy, V.L., G. Gianconne, et al. (2011). "American Society of Clinical Oncology Provisional Clinical Opinion: epidermal growth factor receptor (EGFR) mutation testing for patients with advanced non-small cell lung cancer considering first-line EGFR tyrosine kinase inhibitor therapy." J. Clin. Oncol. 29(15):2121-2127. View Citation Online
26	erlotinib, gefitinib	EGFR	Fukuoka, M., T.S.K. Mok, et al. (2011). "Biomarker analyses and final overall survival results from a phase III, randomized, open-label, first-line study of gefitinib versus carboplatin/paclitaxel in clinically selected patients with advanced non-small-cell lung cancer in Asia (IPASS). J. Clin. Oncol. DOI: 10.1200/JCO.2010.33.4235. View Citation Online
27	dabrafenib and trametinib combination therapy, vemurafenib	BRAF	Hyman, D.H., J. Baselga, et al. (2015). "Vemurafenib in Multiple Nonmelanoma Cancers with BRAF V600 Mutations." NEJM 373(8):726-736. View Citation Online
28	dabrafenib and trametinib combination therapy, vemurafenib	BRAF	Planchard, D., B.E. Johnson, et al. (2016). "An open-label phase II trial of dabrafenib (D) in combination with trametinib (T) in patients (pts) with previously treated BRAF V600E-mutant advanced non-small cell lung cancer (NSCLC; BRF113928)." J Clin Oncol 34: 15_suppl, 107-107. View Citation Online
29	dabrafenib and trametinib combination therapy, vemurafenib	BRAF	Planchard, D., B.E. Johnson, et al. (2017). "Dabrafenib plus trametinib in patients with previously untreated BRAFV600E-mutant metastatic non-small-cell lung cancer: an open-label, phase 2 trial." Lancet Oncol 18(1):1307-1316. View Citation Online
30	dabrafenib and trametinib combination therapy, vemurafenib	BRAF	Subbiah, V., D.M. Hyman, et al. (2019). "Efficacy of Vemurafenib in Patients With Non-Small-Cell Lung Cancer With BRAF V600 Mutation: An Open-Label, Single-Arm Cohort of the Histology-Independent VE-BASKET Study." JCO Precis Oncol 2019: 3, 1-9. View Citation Online
31	pembrolizumab	PD-L1 (22c3)	Herbst, R., P. Baas, et al. (2020). "Long-Term Outcomes and Retreatment Among Patients With Previously Treated, Programmed Death-Ligand 1-Positive, Advanced Non-Small-Cell Lung Cancer in the KEYNOTE-010 Study" J Clin Oncol. 2020 Feb 20;JCO1902446. doi: 10.1200/JCO.19.02446 View Citation Online

References

#	Drug	Biomarker	Reference
32	pembrolizumab	PD-L1 (22c3)	Gadgeel, S., C. Garassino, et al. (2020). "Updated Analysis From KEYNOTE-189: Pembrolizumab or Placebo Plus Pemetrexed and Platinum for Previously Untreated Metastatic Nonsquamous Non-Small-Cell Lung Cancer" J Clin Oncol. 2020 Mar 9;JCO1903136. doi: 10.1200/JCO.19.03136 View Citation Online
33	pembrolizumab	PD-L1 (22c3)	Reck, M., JR Brahmer, et al. (2019). "Updated Analysis of KEYNOTE-024: Pembrolizumab Versus Platinum-Based Chemotherapy for Advanced Non-Small-Cell Lung Cancer With PD-L1 Tumor Proportion Score of 50% or Greater." J Clin Oncol. 37(7):537-546 View Citation Online
34	pembrolizumab	PD-L1 (22c3)	Mok, T.S., KEYNOTE-042 Investigators, et al. (2019). "Pembrolizumab versus chemotherapy for previously untreated, PD-L1-expressing, locally advanced or metastatic non-small-cell lung cancer (KEYNOTE-042): a randomised, open-label, controlled, phase 3 trial." Lancet. 393(10183):1819-1830 View Citation Online